



ECLIPSE BIOINNOVATIONS

ECLIP: PREPARATION AND UV CROSSLINKING OF SUSPENSION CELLS

Required materials:

1. UV crosslinker with 254-nm wavelength UV bulbs (*UVP CL-1000 Ultraviolet Crosslinker or equivalent*)
2. Liquid nitrogen (sufficient to submerge tubes in appropriate container)
3. 1x DPBS (Corning cat #21-031-CV or equivalent)
4. Trypan blue stain (*Thermo Fisher Scientific, Cat# 15250-061 or other equivalent live cell counting assay*)
5. Standard cell counting system (hemocytometer or automated cell counter)

Cell viability validation (prior to crosslinking):

- a. Use Trypan blue stain (Thermo Fisher Scientific, Cat# 15250-061) or other equivalent live cell counting assay to assay cell viability
- b. Cell viability should be > 95% to ensure intact RNA

Preparation of suspension cells:

- a. Pool all cells per biosample (if multiple plates):
- b. Transfer cells with media to 50mL conical tube(s)
- c. Centrifuge at 200g for 5 minutes at room temperature
- d. Aspirate spent media

Wash cells:

- a. Resuspend the pellet(s) in 25 mL of 1x DPBS at room temperature.
- b. Count cell concentration (either with automated cell counter or hemocytometer)
- c. Spin down remaining sample in 50ml conical tube(s) at 200g for 5 minutes at room temperature.
- d. Aspirate supernatant
- e. Resuspend cells to no more than 20×10^6 cells per mL

UV crosslinking:

- a. Aliquot at most 60×10^6 cells (re-suspended in 1x DPBS) in at least 3 mL total volume to a standard 10cm tissue culture grade plate.
 - **Note:** Ensure the cells are evenly dispersed and the plate is fully covered (3 mL should be sufficient volume).
- b. Place the tissue culture plate on leveled ice or a cooling block pre-chilled to 4°C
- c. Place the above (plate plus ice or cooling block) into the UV cross-linker.
 - **Notes: Ensure the plate is leveled.**
 - **Remove tissue culture plate lid** for cross-linking.
- d. Cross-link at 254-nm UV with an energy setting of 400 mJoules/cm²



ECLIPSE BIOINNOVATIONS

- **Note:** this is a setting of **4000** on many cross-linkers which display values in 0.1 mJoules/cm²)
- e. After crosslinking is completed, transfer cells to a 50 mL conical tube.
- f. Wash plate once with 7mL of 1x DPBS and add to the same 50 mL tube.
- g. Count cell concentration (either with automated cell counter or hemocytometer)
 - **Note:** ensure cells are re-suspended well before counting
- h. Centrifuge the 50 ml conical tube at 200g for five minutes at room temperature.
- i. Aspirate and discard supernatant.
- j. Resuspend in the desired amount for flash freezing
 - Typically 20×10⁶ cells per mL.
- k. Transfer desired amount into 1.5 mL Eppendorf Safe-Lock Tubes (or equivalent)
 - Typically 1mL of 20×10⁶ cells per mL.
- l. Spin down at 200g for five minutes at room temperature.
- m. Aspirate the supernatant and freeze by submerging the epi-tubes completely in liquid nitrogen.
- n. After frozen (at least thirty seconds), remove from the liquid nitrogen and store at -80°C.